

ARIANA Sequential Compression Comfort Sleeves

The Ariana Kaveh model is designed to work with the Kendall express 9525 SCD and Kendall 700 series SCD systems.

Compatible with Kendall express knee sleeves .



Deep Vein Thrombosis Sleeve or Anti Embolism Sleeve:

It is indicated for prevention of deep vein thrombosis and pulmonary embolism and in general venous thromboembolism for high risk patients.

Circumferential pressure therapy is achieved by tightening of the cuff around calf and thigh.

It facilitates blood circulation with sequential pressure applied to three chambers from below the knee and thigh. This sequential movement is similar to milking action. Each of the three chambers are inflated with a different pressure which facilitates blood circulation and prevention of deep vein thrombosis.

DVT sleeve of Arianamed with Kaveh model is designed to be compatible with Kendall Expresss 9525 SCD Kendall or Covidien SCD 700.



کاف تسهیل کننده گردش خون (DVT Sleeve) یا جوراب ضد آمبولی:

جهت جلوگیری از ایجاد ترومبوز وریدی عمقی (Deep Vein Thrombosis) و آمبولی ریوی (Pulmonary Embolism) و به طور کلی Venous Thromboembolism برای بیمارانی که ریسک بالایی دارند به کار می رود.

فشار درمانی محیطی با بسته شدن کاف به دور ناحیه ساق و ران پا انجام می شود. امکان گردش و حرکت خون با اعمال فشار به صورت ترتیبی به سه محفظه از پایین ساق پا تا ران مهیا می شود. این حرکت به صورت باد شدن ترتیبی و همانند عمل شیر دوشی انجام می شود. هر کدام از سه محفظه با فشاری متفاوت باد شده و موجب تسهیل گردش خون و جلوگیری از ایجاد لخته در عمق عروق می شود.

کاف تسهیل کننده گردش خون آریانا با مدل کاوه برای کار با دستگاه Kendall Express SCD 9525 یا Covidien سری SCD 700 طراحی شده است.

R.code

DVT_20